

Product Datasheet

LILRA3 antibody (Biotin) (orb418294)

Description	LILRA3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Leukocyte immunoglobulin-like receptor subfamily A member 3 protein (222-313AA)
Target	LILRA3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N6C8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.