

Product Datasheet

FAM84A antibody (FITC) (orb418337)

Description	FAM84A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Protein FAM84A protein (15-144AA)
Target	FAM84A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96KN4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.