

Product Datasheet

CLPTM1L antibody (Biotin) (orb418378)

Description	CLPTM1L antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Cleft lip and palate transmembrane protein 1-like protein (32-172AA)
Target	CLPTM1L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96KA5
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.