



Product Datasheet

PDP1 antibody (orb418511)

Description	PDP1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human [Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 1, mitochondrial protein (1-154AA)
Target	PDP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P0J1
Dilution Range	WB:1:500-5000, IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of SH-SY5Y whole ...