


Product Datasheet

2-hydroxyisobutyryl-HIST1H3A (K122) antibody (orb418596)

Description	2-hydroxyisobutyryl-HIST1H3A (K122) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide sequence around site of 2-hydroxyisobutyryl-Lys (122) derived from Human Histone H3.1
Target	HIST1H3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P68431
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
Jurkat
whole c...