



Product Datasheet

H3F3A (Ab-79) antibody (orb418607)

Description	H3F3A (Ab-79) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, WB
Immunogen	Peptide sequence around site of Lys (79) derived from Human Histone H3.3
Target	H3F3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P84243
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
Hela (4*10...



Western blot
analysis of Hela
whole cel...