

Product Datasheet

Di-methyl-Histone H3.1 (K4) antibody (orb418637)

Description	Rabbit monoclonal antibody to Di-methyl-Histone H3.1 (K4)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, WB
Immunogen	A synthesized peptide
Target	HIST1H3A
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P68431
Dilution Range	WB:1:500-1:2000, IF/ICC:1:50-1:500
Expiration Date	12 months from date of receipt.