

## Product Datasheet

### Mono/Di/Tri-methyl-Histone H3.1 (K14) antibody (orb418648)

|                            |                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit monoclonal antibody to Mono/Di/Tri-methyl-Histone H3.1 (K14)                               |
| <b>Species/Host</b>        | Rabbit                                                                                            |
| <b>Reactivity</b>          | Human                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                      |
| <b>Tested Applications</b> | ELISA                                                                                             |
| <b>Immunogen</b>           | A synthesized peptide                                                                             |
| <b>Target</b>              | HIST1H3A                                                                                          |
| <b>Preservatives</b>       | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Form/Appearance</b>     | Liquid                                                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                     |
| <b>Note</b>                | For research use only                                                                             |
| <b>Isotype</b>             | Rabbit IgG                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                        |
| <b>Uniprot ID</b>          | <b>P68431</b>                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                   |