

Product Datasheet

Hydroxyl-Histone H2A type 1-B/E (Y39) antibody (orb418652)

Description	Rabbit monoclonal antibody to Hydroxyl-Histone H2A type 1-B/E (Y39)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthesized peptide
Target	HIST1H2AB
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P04908
Dilution Range	Recommended dilution: WB:1:500-1:2000
Expiration Date	12 months from date of receipt.