

Product Datasheet

Histone H2B type 1-K antibody (orb418674)

Description	Rabbit monoclonal antibody to Histone H2B type 1-K
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A synthesized peptide
Target	HIST1H2BK
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	O60814
Dilution Range	WB:1:5000-1:10000, IHC-P:1:50-1:500
Expiration Date	12 months from date of receipt.