

Product Datasheet

Human CHGA protein (orb418711)

Description	Recombinant Human Chromogranin-A
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	Tag Info: N-terminal GST-taggedExpression Region: 224-457aaSequence Info: Extracellular Domain
Protein Sequence	QAKREEEEEEEAEAGEEAVPEEEGPTTVLNPSPSLGYKEIRKGESRSEALAVD GAGKPGAEAAQDPEGKGEQEHSQQKEEEEEMAVVPQGLFRGGKSGELEQEEE RLSKEWEDSKRWSKMDQLAKELTAEKRLEGQEEEEEDNRDSSMKLSFRARAYGF RGPGPQLRRGWRPSSREDSLEAGLPLQVRGYPEEKKEEEGSANRRPEDQELESL SAIEAELEKVAHQALRRG
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	53.4 kDa
Uniprot ID	P10645
Expiration Date	6 months from date of receipt.