

Product Datasheet

Human HNRNPA2B1 protein (orb418729)

Description	Recombinant Human Heterogeneous nuclear ribonucleos A2/B1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	Tag Info: N-terminal GST-taggedExpression Region: 1-249aaSequence Info: Cytoplasmic Domain
Protein Sequence	MEREKEQFRKLFIGGLSFETTEESLRNYEQWGKLTDCVVMRDPASKRSRGFGF VTFSSMAEVDAAAMAARPHSIDGRVVEPKRAVAREESGKPGAHVTVKKLFVGGIK EDTEEHHLRDYFEEYGKIDTIEIITDRQSGKKRGFGFVTFDDHDPVDKIVLQKYHTI NGHNAEVRKALSRQEMQEDLEVAILEVAPVMEEEEEDMVVEDLDMATRVGATE VVMTTMEEEIMEVEITMILEIITSNLLTTVQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	55.4 kDa
Uniprot ID	I6L957
Expiration Date	6 months from date of receipt.