

Product Datasheet

Human KPNA1 protein (orb418732)

Description	Recombinant Human Importin subunit alpha-1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-tagged Expression Region: 8-538aa Sequence Info: Partial
Protein Sequence	NFRLKSYKNKSLNPDEMRRRREEEGLQLRKQKREEQLFKRRNVATAEEETEEEV MSDGGFHEAQISNMEMAPGGVITSDMIEMIFSKSPEQQLSATQKFRKLLSKEPNP PIDEVISTPGVVARFVEFLKRKENCTLQFESAWVLTNIASGNSLQTRIVIQAGAVPI FIELLSSEFEDVQEQAVWALGNIAGDSTMCRDYLDCNLPPLLQLFSKQNRRLTM TRNAVWALSNLCRGKSPPEFAKVSPCLNVLSWLLFVSDTDVLADACWALSYS DGPNDKIQAVIDAGVCRRLVELLMHNDYKVVSPALRAVGNIVTGDDIQTQVILNC SALQSLHLHLLSSPKESIKKEACWTISNITAGNRAQIQTVIDANIFPALISILQTAEFRT RKEAAWAITNATSGGSAEQIKYLVELGCIKPLCDLLTVMDSKIVQVALNGLNLR LGEQEAKRNGTGINPYCALIEEAYGLDKIEFLQSHENQEIQKAFDLIEHYFGTED EDSSIAPQVDLNQQQYIFQQCEAPMEGFQL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	63.5 kDa
Uniprot ID	P52294
Expiration Date	6 months from date of receipt.