

Product Datasheet

Plant Cyanovirin-N homolog protein (orb418839)

Description	Recombinant Triangle waterfern Cyanovirin-N homolog
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Tag Info: NO-taggedExpression Region: 28-142aaSequence Info: Full Length
Protein Sequence	QCNFANSCTGVELYGYILRGDCINEDGHPHATSINLNYYIGNDNGRLEYPGESFG SSCVKTALNDGHTLTASCKGADGQYHDSSMDLNYYVVGNSYGYMEPCRASNAD HVLKSSSE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	12.4 kDa
Uniprot ID	P86326
Expiration Date	6 months from date of receipt.