

Product Datasheet

Human DYNC2LI1 protein (orb418869)

Description	Recombinant Human Cytoplasmic dynein 2 light intermediate chain 1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	Tag Info: N-terminal GST-tagged Expression Region: 1-352aa Sequence Info: Full Length
Protein Sequence	MPSETLWEIAKAEVEKRGINGSEGDGAEIAEKVFVFFIGSKNGGKTTIILRCLDRDE PPKPTLALEYTYGRRRAKGHNTPKDIAHFWELGGGTSLLDLISIPITGDTLRTFSLVL VLDLSKPNDLWPTMENLLQATKSHVDKVMKLGKTNAAVSEMQRQKIWNMMPK DHPQDHELIDPFVPLVIIGSKYDVFQDFESEKRKVICKTLRFVAHYYGASLMFTS KSEALLLKIRGVINQLAFGIDKSKSICVDQNKPLFITAGLDSFGQIGSPVPENDIGK LHAHSPMELWKKVYEKLFPPKSINTLKDIKDPARDPQYAENEVDQMRIQKDLELE QYKRSSSKSWKQIELDS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	66.8 kDa
Uniprot ID	Q8TCX1
Expiration Date	6 months from date of receipt.