

Product Datasheet

Mouse DCT protein (orb418893)

Description	Recombinant Mouse L-dopachrome tautomerase
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-SUMO-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 10xHis-SUMO-tagged Expression Region: 24-472aa Sequence Info: Full Length
Protein Sequence	QFPRVCMTLDGVLNKECCPLGPEATNICGFLEGRGQCAEVQTDTRPWSGPYIL RNQDDREQWPRKFFNRTCKCTGNFAGYNCGGCKFGWTGPDCNRKKPAILRRNI HSLTAQEREQFLGALDLAKKSIHPDYVITTQHWLGLLPNGTQPQIANCSVYDFF VWLHYYSVRDTLLGPRPYKAIDFSHQGPAFVTWHRYHLLWLERELQRLTGNES FALPYWNFATGKNECDVCTDELLGAARQDDPTLISRNSRFSTWEIVCDSLDDYN RRVTLCNGTYEGLLRRNKVGRNNEKLPTLKNVQDCLSLQKFDSPFFQNSTFSF RNALEGFDKADGTLDSQVMNLHNLHNSFLNGTNALPHSAANDPVFVVLHSFTD AIFDEWLKRNNPSTDPAWPQELAPIGHNRMYNMVPPFPVTNEELFLTAEQLGYN YAVDLSEEEAPVWSTTLS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	69.7 kDa
Uniprot ID	P29812
Expiration Date	6 months from date of receipt.