

Product Datasheet

E. coli FLIC protein (orb418925)

Description	Recombinant Escherichia coli Flagellin
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-SUMO-tagged Expression Region: 2-498aa Sequence Info: Partial
Protein Sequence	AQVINTNSLSLITQNNINKNQSALSSSIERLSSGLRINSAKDDAAGQAIANRFTSNI KGLTQAARNANDGISVAQTTEGALSEINNQLRVRELTQATTGTNSES DLSSIQ DEIKSRLDEIDRVSGQTQFNGVNV LAKNGSMKIQVGANDNQ TITIDLKQIDAKTL GLDGF SVKNNDTVTTSAPVTA FGATTTNNIKLTGITLSTEAATDTGGTNPASIEGV YTDNGNDYYAKITGGDNDGKYYAVTVANDGTVT MATGATANATVTDANTTKAT TITSGGTPVQIDNTAGSATANLGAVSLVKLQDSKGNDTD TYALKD TNGNLYAAD VNETTGAVSVKTITYTDSSGAASSPTAVKLGGDDGKTEVV DIDGKTYDSADLNG GNLQTGLTAGGEALTAVANGKTTDPLKALDDAIASVDKFRSSLGAVQNR LDSAV TNLNNTTTNLSEAQSRIQDADYATEVSNMSKAQIIQQAGNSVLAKANQVPQQVL SLLQG
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	67.2 kDa
Uniprot ID	P04949
Expiration Date	6 months from date of receipt.