

Product Datasheet

Mouse CES1C protein (orb418954)

Description	Recombinant Mouse Carboxylesterase 1C
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Tag Info: NO-taggedExpression Region: 19-550aaSequence Info: Extracellular Domain
Protein Sequence	HSLLPPVDDTTQGKVLGKYISLEGFEQPVAVFLGVPFAKPLGLSLRFAPPQPAEP WSVFKNATSYPPMCSQDAGWAKILSDMFSTEKEILPLKISEDCLYLNISPADLTK SSQLPVMVWIHGGGLVIGGASPYNGLALSAHENVVVTIQYRLGIWGLFSTGDE HSPGNWAHL DQLAALRWVQDNIANFGGNPDSVTIFGESSGGISVSVLVLSPLGK DLFHRAISESGVVINTNVGKKNIQAVNEIATLSQCNDTSSAAMVQCLRQKTESEL LEISGKLVQYNISLSTMIDGVVLPKAPPEILAEKSFNTVPYIVGFNKQEFGWII PMM LQNLLPEGKMNEETASLLRRFHSSELNISESMIPAVIEQYLRGVDDPAKKSELILD MFGDIFFGIPAVLLSRSLRDAGVSTYMYEFRYRPSFVSDKRPQTVEGDHGEIFF VFGAPLLKEGASEEETNLSKMVMKFWANFARNGNPNNGEGLPHWPEYDEQEGYL QIGATTQQAQRLKAEVAFWTELLAKNPPEPTDPTHEH
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	58.6 kDa
Uniprot ID	P23953
Expiration Date	6 months from date of receipt.