

Product Datasheet

Human p53 protein (orb418963)

Description	Recombinant Human Cellular tumor antigen p53
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-tagged Expression Region: 1-393aa Sequence Info: Full Length
Protein Sequence	MEEPQSDPSVEPPLSQETFSDLWKLLPENNVLSPPLSQAMDDLMLSPDDIEQWF TEDPGPDEAPRMPEAAAPPVAPAPAAPTPAAPAPAPSWPLSSSVPSQKTYQGSYG FRLGFLHSGTAKSVTCTYSPALNKMFCQLAKTCPVQLWVDSTPPPGTRVRAMAI YKQSQHMTEVVRCPHHERCSDSDGLAPPQHILRVEGNLRVEYLDDRNTFRHS VVVPYEPPEVGSDCTTIHYNMCMSSCMGGMNRRPILITLEDSSGNLLGRNSF EVRVCACPGRRDRRTEENLRKKGEPPHELPPGSTKRALPNNTSSSPQPKKKPLD GEYFTLQIRGRERFEMFRELNEALELKDAQAGKEPGGSRAHSSHLKSKKGQSTS RHKKLMFKTEGPDSD
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	45.7 kDa
Uniprot ID	P04637
Expiration Date	6 months from date of receipt.