

Product Datasheet

Bacterial CLFA protein (orb418971)

Description	Recombinant Staphylococcus aureus Clumping factor A
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Tag Info: NO-taggedExpression Region: 228-558aaSequence Info: Full Length
Protein Sequence	GKDITNQLTNVTVGIDSGDTVYPHQAGYVKLNYGFSVPNSAVKGDTFKITVPKEL NLNGVTSTAKVPPIMAGDQVLANGVIDSDGNVIYTFDYYVDTKENVTANITMPAY IDPENVTKTGNVTLTTGIGSTTANKTVLVDYEKYGKFYNLSIKGTIDQIDKTNNNTYR QTIYVNPSPGDNVIAPVLTGNLKPNTDSNALIDQQNTSIKVYKVDNAADLSESYFV NPENFEDVTNSVNITFPNPNQYKVEFNTDQITTPYIVVVNGHIDPNSKGDALRL STLYGYDSRFVWRSMWDNEVAFNNGSGSGDGIDKPVVPEQPDEPGEIEPIPE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	36.1 kDa
Uniprot ID	Q6GB45
Expiration Date	6 months from date of receipt.