

Product Datasheet

Animal Limulus clotting factor C protein (orb418994)

Description	Recombinant Tachypleus tridentatus Limulus clotting factor C,partial
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-taggedExpression Region: 763-1019aaSequence Info: Full Length
Protein Sequence	IWNGNSTEIGQWPWQAGISRWLADHNMWFLQCGGSLLNEKWIVTAAHCVTYS ATAEIIDPSQFKIYLGKYRDDSRRDDYVQVREALEIHVNPNYDPGNLNFIDIALIQL KTPVTLTTRVQPICLPTDITTREHLKEGTLAVVTGWGLNENNTYSEMIQQAVLPV VAASTCEEYKEADLPLVTENMFCAGYKKGRYDACSGDSGGPLVFADDSRTE RRWVLEGIVSWGSPSGCGKANQYGGFTKVNVLFSWIRQFI
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Yeast
MW	32.3 kDa
Uniprot ID	P28175
Expiration Date	6 months from date of receipt.