

Product Datasheet

Mouse TLR4 protein (orb419009)

Description	Recombinant Mouse Toll-like receptor 4
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-tagged Expression Region: 26-638aa Sequence Info: Full Length
Protein Sequence	NPCIEVVPNITYQCMDQKLSKVPDDIPSSTKNIDLSFNPLKILKSYSFSNFSELQWL DLSRCEIETIEDKAWHGLHLSNLILTGNPIQSFSFGSFGSLTLENLVAVETKLAS LESFPIGQLITLKKLNVAHNFIHSCKLPAYFSNLTNLVHVDLSYNYIQTITVNDLQFL RENPQVNLSDMSLNPIDFIQDQAFQGIKHELTLRGNFNSSNIMKTCLQNLAGL HVHRLILGEFKDERNLEIFEPSIMEGLCDVTIDEFRLTYTNDFSDDIVKFHCLANVS AMSLAGVSIKYLEDPKHFKWQSLSIIRCQLKQFPTLDLPFLKSLTLTMNKGSIKFK KVALPSLSYLDLSRNALSFSGCCSYSDLGTNSLRHLDLSFNNGAIIIMSANFMGLEEL QHLDFQHSTLKRVTESAFSLSEKLLYLDISYTNNTKIDFDGIFLGLTSLNTLKMAGN SFDNTLSNVFANTTNLTFLDLSKCQLEQISWGVFDTLHRLQLLNMSHNNLLFLD SSHYNQLYSLTDCSFNRIETSKGILQHFPKSLAFFNLNNSVACICEHQKFLQW VKEQKQFLVNVEQMTCATPVEMNTSLVLDNFNNSTCYMYKTIISVSVVS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	71.5 kDa
Uniprot ID	Q9QUK6
Expiration Date	6 months from date of receipt.