

Product Datasheet

Mouse FLII protein (orb419023)

Description	Recombinant Mouse flightless-1 homolog
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-B2M-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-B2M-tagged Expression Region: 495-827aa Sequence Info: Partial
Protein Sequence	VGQLPGLTIWQIENFVPLVEEAFHGKFYEADCYIVLKTFLDDSGSLNWEIYYWIG GEATLDKKACSAIHAVNLRNYLGAECRTVREEMGDESEEFLLQVFDNDISYIEGGT ASGFYTVEDTHYVTRMYRVYGKKNIKLEPVPLKGSSLDPRFVFLDQGLDIYVWR GAQATLSNTTKARLFAEKINKNERKGKAEITLLVQGGQEPGFWDVLGGEPSEIKN HVPDDFWPPQPKLYKVGGLGYLELPQINYKLSVEHKKRPKVELMPGMRLQLSLL DTRCVYILDCWSDVFIWLGRKSPRLVRAAALKGQELCGMLHRPRHTTVVSRSL GTE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	52 kDa
Uniprot ID	Q9JJ28
Expiration Date	6 months from date of receipt.