

Product Datasheet

Xenopus GATA1-A protein (orb419027)

Description	Recombinant Xenopus laevis GATA-binding factor 1-A
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged Expression Region: 1-359aa Sequence Info: Full Length of Mature Protein
Protein Sequence	MDYTTLTQDPDPNYTESGLASTSEDSQFLYGLGGESSPGHYGGAVSSRAVGG FRHSPVFQTFPLHWPETSAGIPSNLTAYGRSTGTLSFYPSAASALGPITSPPLYSAS SFLGSAAPPAEREGSPKFLETLKTERASPLTSDLLPLEPRSPSILQVGYIGGGGQEF SLFQSTEDRECVNCGATVTPLWRRDMSGHYLCNACGLYHKMNGQNRPLIRPKK RLIVSKRAGTQCSNCHTSTTTLWRRNASGDPVCNACGLYYKLHNVNRPLTMKKE GIQTRNRKVSSRSKKKKQLDNPFEPPKAGVEEPSYPFGPLL FHGQMPPMGHMI NPPHHFLQSPRISHSAPAVSYRQAASGVTPP
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	59 kDa
Uniprot ID	P23767
Expiration Date	6 months from date of receipt.