

Product Datasheet

Bacterial PAC protein (orb419097)

Description	Recombinant Streptomyces alboniger Puromycin N-acetyltransferase
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-taggedExpression Region: 1-199aaSequence Info: Full Length
Protein Sequence	MTEYKPTVRLATRDDVPRAVRTLAAAFADYPATRH TVDPDRHIERVTELQELFLT RVGLDIGKVVVADDGAAVAVWTTPE SVEAGAVFAEIGPRMAELSGSRLAAQQQ MEGLLAPHRPKPAWFLATVGVSPDHQKG LGS AVVLPGV EAAERAGVPAFLE TSAPRNLPFYERLGFTVTADVEVPEGPRTWCMTRKPGA
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	25.5 kDa
Uniprot ID	P13249
Expiration Date	6 months from date of receipt.