

Product Datasheet

Bovine SV2A protein (orb419102)

Description	Recombinant Bovine Synaptic vesicle glyco 2A
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-tagged and C-terminal Myc-tagged Expression Region: 1-169aa Sequence Info: Full Length
Protein Sequence	MEEGFRDRAAFIRGAKDIAKEVKKHATKKVVKGLDRVQDEYSRRSYSRFEEEDD DDDFPAPADGYRGEQAQDEEEGGASSDATEGHDEDEIYEGEYQGIPRAESG GKGERMADGAPLAGVRGGLGDGEGPPGGRGEAQRRKEREELAQQYEAILREC GHGRFQWTLY
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	22.3 kDa
Uniprot ID	Q29397
Expiration Date	6 months from date of receipt.