

Product Datasheet

Fungi PHYA protein (orb419137)

Description	Recombinant Aspergillus niger 3-phytase A
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-taggedExpression Region: 24-467aaSequence Info: Full Length of Mature Protein
Protein Sequence	ASRNQSSCDTVDQGYQCFSETSHLWGQYAPFFSLANESVISPEVPAGCRVTFAQ VLSRHGARYPTDSKGKKYSALIEEIQQNATTFDGKYAFLKTYNYSLGADDLTPFG EQELVNSGIKFYQRYESLTRNIVPFIRSSGSSRVIASGKKFIEGFQSTKLKDPRAQP GQSSPKIDVVISSESSNNTLDPGTCTVFEDSELADTVEANFTATFVPSIRQRLEN DLSGVTLTDEVTYLMDMCSFDTISTSTVDTKLSPFCDLFTHDEWINYDYLQSLK KYYGHGAGNPLGPTQGVGYANELIARLTHSPVHDDTSSNHTLDSSPATFPLNST LYADFSHDNGIISILFALGLYNGTKPLSTTTVENITQTDGFSSAWTVPFASRLYE MMQCQAEQEPLVRVLVNDRVVPLHGCPVDALGRCTRDSFVRGLSFARSGGDW AECFA
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	52.8 kDa
Uniprot ID	P34752
Expiration Date	6 months from date of receipt.