

Product Datasheet

Human FBXL2 protein (orb419149)

Description	Recombinant Human F-box/LRR-repeat 2
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 6xHis-taggedExpression Region: 1-423aaSequence Info: Partial
Protein Sequence	MVFSNNDEGLINKKLPKELLRLIFSLDIVTLRCQAQISKAWNIALDGSNWQRID LFNFQTDVEGRVVENISKRCGGFLRKLSLRGCIGVGDSSSLKTFAQNCRNIEHLNL NGCTKITDSTCYLSRFCSKLKHLDTSCVSITNSSLKGISEGCRNLEYLNLSWCD QITKDIEALVRGCRGLKALLRGCTQLEDEALKHIQNYCHELVSLNLQSCSRITD EGVVQICRGCHRLQALCLSGCSNLTASLTALGLNCPRLQILEAARCSHLTDAGF TLLARNCHELEKMDLEECILITDSTLIQLSIHCPKLQALSLSHCELITDDGILHLSNS TCGHERLRVLELDNCLITDVALEHLENCRGLERLELYDCQVTRAGIKRMRAQL PHVKVHAYFAPVTPPTAVAGSGQRLCRCCVIL
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Yeast
MW	49.1 kDa
Uniprot ID	Q9UKC9
Expiration Date	6 months from date of receipt.