

Product Datasheet

Mouse ISM1 protein (orb419157)

Description	Recombinant Mouse Isthmin-1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 6xHis-tagged
Note	For research use only
Application notes	Tag Info: C-terminal 6xHis-tagged Expression Region: 30-454aa Sequence Info: Full Length of Mature Protein
Protein Sequence	SDRQDAAAGNNLNLESdstsetSFPLSKEAPEEHQVVHQPFPRQRFPPEtGHP SLQRDGP RSFLLDLPNFPDLSKADINGQNPNIQVTIEVVDGPDSEAEKDQHPENK PSWSLPAPDWRAWWQRSLSLARTNSGDQDDKYDSTSDDSNFLSVPRGWDRP APGHRTFETKEQPEYDSTDGEGDWLSVCSVTCGNGNQKRTRSCGYACIATE SRTC DRPNCPGIEDTFR TAATEVSLLAGSEEFNATKLFEVDMDS CERWMSCKSE FLKKYMHKVINDLPSCPCSYPTVAYSTADIFDRIKRKDFRWKDASGPKEKLEIYK PTARYCIRSM LSLESTTLAAQHCCYGD NMQLITRGKGAGTPNLISTEFSAELHYKV DVL PWIICKGDWSRYNEARPPNNGQKCTESPSDEDYIKQFQEAREY
Purity	Greater than 85% as determined by SDS-PAGE.
Source	Yeast
MW	50.1 kDa
Uniprot ID	A2ATD1
Expiration Date	6 months from date of receipt.