

Product Datasheet

Fungi 6GAL protein (orb419166)

Description	Recombinant Hypocrea rufa Endo-beta-1,6-galactanase
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 10xHis-tagged Expression Region: 21-479aa Sequence Info: Full Length
Protein Sequence	DTTLTIDPTSNWGTWEGWGVSLAWWAKAFGNRDDLASVFFSRNNQAVNGQT LPGLGFNIVRYNAGACSNNSYDGSTMVVSPNIKPSRQMDGFWLDWASSDPSSS SWNWNVDANQRAMLQKAKANGANIFELFSNPMWWMCNNHNPSGSGSSDN LQSWNYQNHAVYLADIAQHAQQSWRIQFSVEAFNEPSSSWWTAEGTQEGCH FDVSTMATVIGYLNTELSSRGLSSFVASSDENTYDLAISTWQGFNSSTRNIVKRIN VHGYQDGGGRRDLYSLASQAGKRLWNSEYGDSDASGKSMYQNLLEDFTWLH PTAWVYWQAIDGAGWGLIVGDNDNLTLSSASTKYFVLAQLTRHIRQGMQILTP DVNTAVAYDAGSQKLIVTANWGSAQTTITFDLTRARTAGSNGATVPRWSTQTG GGDQYRSYTDTKINNGKFSASFSSGQVQTFEVSGVVLQ
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	56.1 kDa
Uniprot ID	Q76FP5
Expiration Date	6 months from date of receipt.