

Product Datasheet

Primate IL5 protein (orb419312)

Description	Recombinant Rhesus Macaque Interleukin-5 active
Endotoxins	Less than 1.0 EU/μg as determined by LAL method.
Preservatives	Lyophilized from a 0.2 μm filtered PBS, pH 7.4
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Tag Info: NO-taggedExpression Region: 20-134aaSequence Info: Partial
Protein Sequence	IPTEIPASALVKETLALLSTHRTLLIANETLRIPVPVHKNHQLCTEEIFQGIGTLESQT VQGGTVERLKFKNLSLIKKYIGGQKKKCGEERRRVNQFLDYLQEFLGVMNTEWIIIE S
Purity	> 98% as determined by SDS-PAGE and HPLC.
Source	E.Coli
MW	13.1 kDa
Uniprot ID	P48093
Expiration Date	6 months from date of receipt.

