

Product Datasheet

Viral Glycoprotein protein (orb419329)

Description	Recombinant Rabies virus Glyco G
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Tag Info: N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged Expression Region: 20-459aa Sequence Info: Full Length
Protein Sequence	KFPIYTIPDKLGPWSPIDIHHLSCPNLVEDEGCTNLSEFSYMEKLVGYISAIKVN GFTCTGVVTEAETYNFVGIVTTTFKRKHFRPTPDACRAAYNWKMAGDPRYEES LHNYPDPYHWLRTVRTTKESLIISPSVTDLPYDKSLHSRVFPGGKCSGITVSSTY CSTNHDYTIWMPENPRPRTPCDIFTNSRGKRASKGNKTCGFVDERGLYKSLKGA CRLKLCGVLGLRLMDGTWVAMQTSDETKWCPPDQLVNLHDFRSDEIEHLVVEE LVKKREECLDALESIMTTKSVSFRRLSHLRKLVPGFGKAYTIFNKTLMEADAHYKS VRTWNEIIPSKGCLKVGGRCRPHVNGVFFNGIILGPDGHVLIPEMQSSLLQQHME LLKSSVIPLMHPLADPSTVFKEGDEAEDFVEVHLPDVYKQISGVDLGLPNWGKY
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	69.7 kDa
Uniprot ID	O92284
Expiration Date	6 months from date of receipt.