


Product Datasheet

alpha-Linolenic Acid (orb419444)

Description	alpha-Linolenic Acid
Storage	Storage Temp: -20°C
Note	For research use only
Purity	≥98%
MW	278.44
CAS Number	463-40-1
SMILES	<chem>CC/C=C\C/C=C\C/C=C\CCCCCCCC(=O)O</chem>
Solubility (25°C)	Soluble in methanol, ethanol, ethyl acetate. Slightly soluble in chloroform.
Hazard Information	UN Number: Not dangerous goods
Expiration Date	12 months from date of receipt.



Chemical structure of alpha-Linolenic Ac...