

Product Datasheet

Human IL-2 protein (orb419913)

Description	Recombinant of human IL-2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 10mM Acetata-Na, 5% Trehaiose, pH 4.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 10mM Acetata-Na, 5% Trehaiose, pH 4.5 .
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Application notes	C145S
Protein Sequence	Pro22-Thr153(Cys145Ser)
Source	E.coli
Uniprot ID	P60568
Expiration Date	6 months from date of receipt.