

Product Datasheet

Human MMP8 protein (orb419930)

Description	Recombinant of human MMP8 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 50mM Tris-HCl,10mM CaCl ₂ ,150mM NaCl,pH7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mM Tris-HCl,10mM CaCl ₂ ,150mM NaCl,pH7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	C-10His
Protein Sequence	Phe21-Gly467
Source	Human Cells
Uniprot ID	P22894
Expiration Date	6 months from date of receipt.