

Product Datasheet

Rat CD73 protein (orb419936)

Description	Recombinant of rat CD73 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 120mM NaCl, 4mM CaCl ₂ , 20% Glycerol, pH 7.5.
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	C-6His
Protein Sequence	Trp29-Lys549
Source	Human Cells
Uniprot ID	P21588
Expiration Date	6 months from date of receipt.