

Product Datasheet

Mouse CD262 protein (orb419953)

Description	Recombinant of mouse CD262 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-AVI&His
Note	For research use only
Application notes	C-Avi-6His
Protein Sequence	Asn53-Ser177
Source	Human Cells
Uniprot ID	Q9QZM4
Expiration Date	6 months from date of receipt.