

Product Datasheet

Mouse LR3-IGF1 protein (orb419958)

Description	Recombinant of mouse LR3-IGF1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mMTris,150mMNaCl,1mMEDTA,pH8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mMNaAC,ph4.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	C-6His
Protein Sequence	Gly49-Ala118
Source	E.coli
Uniprot ID	P05017
Expiration Date	6 months from date of receipt.