

Product Datasheet

Human Sting protein (orb419962)

Description	Recombinant of human Sting protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His&SUMO
Note	For research use only
Application notes	N-Sumo-6His
Protein Sequence	Val155-Val341
Source	E.coli
Uniprot ID	Q86WV6
Expiration Date	6 months from date of receipt.