

Product Datasheet

Human Troponin I type 3 protein [Out of stock] (orb419963)

Description	Recombinant of human Troponin I type 3 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, 2mM EDTA, 5% Trehalose, pH 7.4.
Tag	N-His
Note	For research use only
Application notes	N-6His
Protein Sequence	Met1-Ser210
Uniprot ID	Q6FGX2
Expiration Date	6 months from date of receipt.