

Product Datasheet

Human INDOL1 protein (orb419964)

Description	Recombinant of human INDOL1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 10%Glycerol, 1mM EDTA, 250mM NaCl,pH8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 10% Glycerol,1mMEDTA,250mM NaCl,pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	C-6His
Protein Sequence	Met14-Gly420
Source	E.coli
Uniprot ID	Q6ZQW0
Expiration Date	6 months from date of receipt.