

Product Datasheet

Human NNE protein (orb419966)

Description	Recombinant of human NNE protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 250mM NaCl, 1mM MgSO4,20%glycerol□Ph7.5 .
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 250mM NaCl, 1mM MgSO4,20% Glycerol,pH 7.5 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	C-6His
Protein Sequence	Met1-Lys434
Source	E.coli
Uniprot ID	P06733
Expiration Date	6 months from date of receipt.