

Product Datasheet

Human TNNC1 protein (orb419967)

Description	Recombinant of human TNNC1 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,100mM NaCl,1mM DTT,10% glycerol,pH8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris,100mM NaCl,1mM DTT,10% Glycerol,pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	N-6His
Protein Sequence	Met1-Glu161
Source	E.coli
Uniprot ID	P63316
Expiration Date	6 months from date of receipt.