

Product Datasheet

Human IL-2 protein (orb419969)

Description	Recombinant of human IL-2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	Mutant
Protein Sequence	Ala21-Thr153(L100F,R101D,L105V,I106V,I112F)
Source	Human Cells
Uniprot ID	P60568
Expiration Date	6 months from date of receipt.