

Product Datasheet

Human sCD4 protein (orb419970)

Description	Recombinant of human sCD4 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, 1mM EDTA, PH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, 1mM EDTA, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Application notes	C-Fc
Protein Sequence	Lys26-Trp390
Source	Human Cells
Uniprot ID	P01730
Expiration Date	6 months from date of receipt.