

Product Datasheet

Human Hemopexin protein (orb419974)

Description	Recombinant of human Hemopexin protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM MES,150mM NaCl,pH5.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM MES, 150mM NaCl, pH 5.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	C-6His
Protein Sequence	Thr24-His462
Source	Human Cells
Uniprot ID	P02790
Expiration Date	6 months from date of receipt.