

Product Datasheet

Human TROP-2 protein (orb419980)

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Recombinant of human TROP-2 protein |
| Endotoxins | Less than 0.1 ng/ug (1 EU/ug) |
| Preservatives | Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4. |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | C-Fc |
| Note | For research use only |
| Application notes | C-Fc |
| Protein Sequence | Thr88-Thr274 |
| Source | Human Cells |
| Uniprot ID | P09758 |
| Expiration Date | 6 months from date of receipt. |