

Product Datasheet

Mouse Serpin E2 protein (orb419982)

Description	Recombinant of mouse Serpin E2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris,150 mM NaCl,pH8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150 mM NaCl, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	C-10His
Protein Sequence	Ser20-Pro397
Source	Human Cells
Uniprot ID	Q07235
Expiration Date	6 months from date of receipt.