

Product Datasheet

Mouse CD39 protein (orb419986)

Description	Recombinant of mouse CD39 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 500mM NaCl, 10% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	C-6His
Protein Sequence	Thr38-Ile478
Source	Human Cells
Uniprot ID	P55772
Expiration Date	6 months from date of receipt.