

Product Datasheet

Mouse CD40LG protein (orb420003)

Description	Recombinant of mouse CD40LG protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 200mM NaCl, 0.1mM EDTA, pH 7.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Application notes	N-6His
Protein Sequence	Met112-Leu260
Source	Human Cells
Uniprot ID	P27548
Expiration Date	6 months from date of receipt.